

USP21 antibody



Catalog Number: 116599

Product name

USP21 antibody

Specificity

Human; other species not tested.

Antibody description

USP21 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, recombinant protein. Positive IHC detected in human heart tissue, human pancreas tissue. Observed molecular weight by Western-blot: 63 kDa

Preparation

This antibody was obtained by immunization of USP21 recombinant protein (Accession Number: NM_012475). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116599(USP21 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116599(USP21 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116599(USP21 Antibody) at dilution of 1:50

USP21 antibody



Catalog Number: 116599

(under 40x lens)