

# USP21 antibody



Catalog Number: 116599

## Product name

USP21 antibody

## Specificity

Human; other species not tested.

## Antibody description

USP21 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, recombinant protein. Positive IHC detected in human heart tissue, human pancreas tissue. Observed molecular weight by Western-blot: 63 kDa

## Preparation

This antibody was obtained by immunization of USP21 recombinant protein (Accession Number: NM\_012475). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116599(USP21 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116599(USP21 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:116599(USP21 Antibody) at dilution of 1:50

# USP21 antibody



Catalog Number: 116599

---

(under 40x lens)