

USP25 antibody



Catalog Number: 116602

Product name

Polyclonal

USP25 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

USP25 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, human cerebellum tissue. Observed molecular weight by Western-blot: 125 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of USP25 recombinant protein (Accession Number: BC015930). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116602(USP25 antibody) at dilution of 1:400