

USP36 antibody



Catalog Number: 116610

Product name

USP36 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

USP36 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HepG2 cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 100-123 kDa

Preparation

This antibody was obtained by immunization of USP36 recombinant protein (Accession Number: BC038983). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116610(USP36 antibody) at dilution of 1:200



IP Result of anti-USP36 (IP:Catalog No:116610, 4ug; Detection:Catalog No:116610 1:200) with Jurkat cells lysate 1200ug.