

# USP36 antibody

Catalog Number: 116610

## Product name

USP36 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

USP36 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HepG2 cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 100-123 kDa

## Preparation

This antibody was obtained by immunization of USP36 recombinant protein (Accession Number: BC038983). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116610(USP36 antibody) at dilution of 1:200



IP Result of anti-USP36 (IP:Catalog No:116610, 4ug; Detection:Catalog No:116610 1:200) with Jurkat cells lysate 1200ug.