

# UNC5CL antibody

Catalog Number: 116617

## Product name

UNC5CL antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

UNC5CL Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of UNC5CL recombinant protein (Accession Number: NM\_173561). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116617(UNC5CL Antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116617(UNC5CL Antibody) at dilution of 1:1000