

UQCR antibody

Catalog Number: 116633

Product name

UQCR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UQCR Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue.

Preparation

This antibody was obtained by immunization of UQCR recombinant protein (Accession Number: NM_006830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:116633(UQCR antibody) at dilution of 1:50 (under 10x lens) 

Immunohistochemical of paraffin-embedded human heart using Catalog No:116633(UQCR antibody) at dilution of 1:50 (under 40x lens)