

# URAT1 antibody



Catalog Number: 116640

## Product name

URAT1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

URAT1 Rabbit Polyclonal antibody. Positive WB detected in PC-13 cells, HEK-293 cells, HepG2 cells, human kidney tissue. Observed molecular weight by Western-blot: 65-70 kDa or 40 kDa

## Preparation

This antibody was obtained by immunization of URAT1 recombinant protein (Accession Number: NM\_144585). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:116640(URAT1 antibody) at dilution of 1:600