

# UBP1 antibody



Catalog Number: 116657

## Product name

UBP1 antibody

## Immunogen

[Human UBP1 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UBP1 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue. Observed molecular weight by Western-blot: 68kd

## Preparation

This antibody was obtained by immunization of UBP1 recombinant protein (Accession Number: NM\_014517). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116657(UBP1 antibody) at dilution of 1:1000