

UBTD2 antibody



Catalog Number: 116665

Product name

UBTD2 antibody

Immunogen

[Human UBTD2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBTD2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of UBTD2 recombinant protein (Accession Number: XM_017010023). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116665(UBTD2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116665(UBTD2 antibody) at dilution of 1:100 (under 10x lens)