

# UBXN1 antibody



Catalog Number: 116668

## Product name

UBXN1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBXN1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HepG2 cells, human brain tissue, human heart tissue, mouse cerebellum tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 43-45kd

## Preparation

This antibody was obtained by immunization of UBXN1 recombinant protein (Accession Number: NM\_001286077). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116668(UBXN1 antibody) at dilution of 1:500