

UBXN11 antibody



Catalog Number: 116670

Product name

UBXN11 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBXN11 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human testis tissue. Positive IP detected in mouse testis tissue. Positive WB detected in PC-3 cells, mouse testis tissue. Observed molecular weight by Western-blot: 48-54kd

Preparation

This antibody was obtained by immunization of UBXN11 recombinant protein (Accession Number: BC038106). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116670(UBXN11 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116670(UBXN11 antibody) at dilution of 1:100 (under 10x lens)

IP Result of anti-UBXN11 (IP:Catalog No:116670, 4ug; Detection:Catalog No:116670 1:500) with mouse testis tissue lysate 4000ug.