

UCMA antibody



Catalog Number: 116678

Product name	Polyclonal
UCMA antibody	Ig Type
Specificity	Rabbit IgG
Human; other species not tested.	Applications
Antibody description	ELISA, WB
UCMA Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.	Dilutions
Preparation	Recommended Dilution:
This antibody was obtained by immunization of UCMA recombinant protein (Accession Number: NM_145314). Purification method: Antigen affinity purified.	WB: 1:500-1:5000
Formulation	Validations
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.	
Storage	Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116678(UCMA Antibody) at dilution of 1:1000
Clonality	