

# USP5 antibody

Catalog Number: 116688

## Product name

USP5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

USP5 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse lung tissue, HeLa cells, human brain tissue, MDA-MB-453s cells, rat brain tissue, rat lung tissue. Observed molecular weight by Western-blot: ~95 kDa

## Preparation

This antibody was obtained by immunization of USP5 recombinant protein (Accession Number: NM\_001098536). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116688(USP5 antibody) at dilution of 1:400



IP Result of anti-USP5 (IP:Catalog No:116688, 4ug; Detection:Catalog No:116688 1:1000) with mouse brain tissue lysate 4000ug.