

UTP14A antibody



Catalog Number: 116694

Product name

UTP14A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UTP14A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, PC-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 88kd

Preparation

This antibody was obtained by immunization of UTP14A recombinant protein (Accession Number: NM_006649). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using UTP14A antibody Catalog No:116694 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116694(UTP14A Antibody) at dilution of 1:1000