

UTP15 antibody



Catalog Number: 116695

Product name

Polyclonal

UTP15 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

UTP15 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 58kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of UTP15 recombinant protein (Accession Number: BC013064). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116695(UTP15 antibody) at dilution of 1:300

