

UXT antibody

Catalog Number: 116700

Product name

UXT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UXT Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 18, 45 kDa

Preparation

This antibody was obtained by immunization of UXT recombinant protein (Accession Number: NM_004182). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116700(UXT antibody) at dilution of 1:400



IP Result of anti-UXT (IP:Catalog No:116700, 3ug; Detection:Catalog No:116700 1:300) with mouse brain tissue lysate 4000ug.