

# VANGL1 antibody



Catalog Number: 116711

---

## Product name

VANGL1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

VANGL1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 60-65 kDa

## Preparation

This antibody was obtained by immunization of VANGL1 recombinant protein (Accession Number: NM\_138959). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116711(VANGL1 antibody) at dilution of 1:300