

# VBP1 antibody



Catalog Number: 116724

## Product name

VBP1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VBP1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IF detected in HepG2 cells. Positive WB detected in rat brain tissue, C6 cells, HepG2 cells, K-562 cells, mouse liver tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 23-25 kDa

## Preparation

This antibody was obtained by immunization of VBP1 recombinant protein (Accession Number: NM\_003372). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, IF, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116724(VBP1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116724(VBP1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using VBP1 antibody Catalog No:116724 at 1:50 dilution

# VBP1 antibody



Catalog Number: 116724

---

and Rhodamine-labeled goat anti-rabbit IgG (red).  
Blue pseudocolor = DAPI (fluorescent DNA dye).

rat brain tissue were subjected to SDS PAGE  
followed by western blot with Catalog  
No:116724(VBP1 Antibody) at dilution of 1:1000

