

VPS11 antibody



Catalog Number: 116770

Product name

VPS11 antibody

Immunogen

[Human VPS11 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS11 Rabbit Polyclonal antibody. Positive FC detected in HepG2 cells. Positive IF detected in HepG2 cells. Positive WB detected in human brain tissue, HEK-293 cells, human heart tissue, human kidney tissue, mouse pancreas tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 108-112 kDa

Preparation

This antibody was obtained by immunization of VPS11 recombinant protein (Accession Number: NM_021729). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, FC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:200-1:2000

IF: 1:10-1:100

Validations





human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116770(VPS11 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using VPS11 antibody Catalog No:116770 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent

VPS11 antibody

Catalog Number: 116770



DNA dye).

IP Result of anti-VPS11 (IP:Catalog No:116770, 4ug; Detection:Catalog No:116770 1:500) with HEK-293 cells lysate 1480ug.

1X10⁶ HepG2 cells were stained with 0.2ug VPS11 antibody (Catalog No:116770, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.