

# VPS37A antibody



Catalog Number: 116785

## Product name

VPS37A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VPS37A Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HeLa cells, MCF-7 cells, mouse liver tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of VPS37A recombinant protein (Accession Number: XM\_017013026). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116785(VPS37A antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116785(VPS37A antibody) at dilution of 1:50 (under 10x lens)