

# VSTM2L antibody

Catalog Number: 116802

## Product name

VSTM2L antibody

## Specificity

Human; other species not tested.

## Antibody description

VSTM2L Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of VSTM2L recombinant protein (Accession Number: NM\_080607). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116802(VSTM2L Antibody) at dilution of 1:1000