

# VSX1 antibody

Catalog Number: 116803

## Product name

VSX1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

VSX1 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, Raji cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of VSX1 recombinant protein (Accession Number: NM\_014588). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:116803(VSX1 Antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116803(VSX1 Antibody) at dilution of 1:5000