

VSX1 antibody



Catalog Number: 116803

Product name

VSX1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

VSX1 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, Raji cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of VSX1 recombinant protein (Accession Number: NM_014588). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:116803(VSX1 Antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116803(VSX1 Antibody) at dilution of 1:5000