

WDR6 antibody



Catalog Number: 116813

Product name

WDR6 antibody

Immunogen

Specificity

human, mouse, rat; other species not tested.

Antibody description

WDR6 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 122kd

Preparation

This antibody was obtained by immunization of WDR6 recombinant protein (Accession Number: BC002826). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116813(WDR6 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using WDR6 antibody Catalog No:116813 at 1:25 dilution and FITC- labeled goat anti-Rabbit IgG (green).