

# WHSC1 antibody

Catalog Number: 116826

## Product name

WHSC1 antibody

## Specificity

Human; other species not tested.

## Antibody description

WHSC1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 70, 150 kDa

## Preparation

This antibody was obtained by immunization of WHSC1 recombinant protein (Accession Number: BC152412). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-WHSC1 (IP:Catalog No:116826, 3ug; Detection:Catalog No:116826 1:300) with HeLa cells lysate 2000ug.