

WHSC1 antibody



Catalog Number: 116826

Product name

WHSC1 antibody

Specificity

Human; other species not tested.

Antibody description

WHSC1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 70, 150 kDa

Preparation

This antibody was obtained by immunization of WHSC1 recombinant protein (Accession Number: BC152412). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-WHSC1 (IP:Catalog No:116826, 3ug; Detection:Catalog No:116826 1:300) with HeLa cells lysate 2000ug.