

# WNT3 antibody

Catalog Number: 116835

## Product name

WNT3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

WNT3 Rabbit Polyclonal antibody. Positive WB detected in SGC-7901 cells, NIH/3T3 cells, Raji cells, SH-SY5Y cells. Positive IP detected in MCF-7 cells. Positive FC detected in HeLa cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of WNT3 recombinant protein (Accession Number: NM\_030753). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, FC, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:116835(WNT3 antibody) at dilution of 1:600



IP Result of anti-WNT3 (IP:Catalog No:116835, 4ug; Detection:Catalog No:116835 1:500) with MCF-7 cells lysate 3200ug.



1X10<sup>6</sup> HeLa cells were stained with 0.2ug WNT3 antibody (Catalog No:116835, red) and control antibody (blue). Fixed with 4% PFA blocked with

# WNT3 antibody



Catalog Number: 116835

---

3% BSA (30 min). Alexa Fluor 488-conjugated

AffiniPure Goat Anti-Rabbit IgG(H+L) with  
dilution 1:1500.