

WAPL; WAPAL antibody



Catalog Number: 116840

Product name

WAPL; WAPAL antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

WAPL; WAPAL Rabbit Polyclonal antibody.
Positive WB detected in HL-60 cells, Jurkat cells.
Observed molecular weight by Western-blot:
140kd

Preparation

This antibody was obtained by immunization of WAPL; WAPAL recombinant protein (Accession Number: BC017393). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116840(WAPL; WAPAL antibody) at dilution of 1:300