

WBSCR16 antibody



Catalog Number: 116850

Product name

WBSCR16 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

WBSCR16 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of WBSCR16 recombinant protein (Accession Number: NM_148842). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116850(WBSCR16 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116850(WBSCR16 antibody) at dilution of 1:50



(under 40x lens)

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116850(WBSCR16 Antibody) at dilution of

WBSCR16 antibody



Catalog Number: 116850

1:600