

# WDFY3 antibody

Catalog Number: 116852

## Product name

WDFY3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

WDFY3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_014991). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116852(WDFY3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116852(WDFY3 antibody) at dilution of 1:50 (under 40x lens)