

WDR4 antibody



Catalog Number: 116866

Product name

WDR4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

WDR4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, transfected cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of WDR4 recombinant protein (Accession Number: NM_001260476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116866(WDR4 antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116866(WDR4 Antibody) at dilution of 1:1000