

# WDR4 antibody

Catalog Number: 116866

## Product name

WDR4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

WDR4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, transfected cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30kd

## Preparation

This antibody was obtained by immunization of WDR4 recombinant protein (Accession Number: NM\_001260476). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116866(WDR4 antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116866(WDR4 Antibody) at dilution of 1:1000