

WDR46 antibody

Catalog Number: 116869

Product name

WDR46 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

WDR46 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HeLa cells. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of WDR46 recombinant protein (Accession Number: NM_005452). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116869(WDR46 antibody) at dilution of 1:400



IP Result of anti-WDR46 (IP:Catalog No:116869, 4ug; Detection:Catalog No:116869 1:400) with Jurkat cells lysate 3200ug.