

YWHAZ antibody



Catalog Number: 116899

Product name

Polyclonal

YWHAZ antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

YWHAZ Rabbit Polyclonal antibody. Positive WB detected in Raji cells, Jurkat cells, mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 28-30 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of YWHAZ recombinant protein (Accession Number: XM_005251061). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116899(YWHAZ Antibody) at dilution of 1:1000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality