

# ZBTB33 antibody



Catalog Number: 116913

## Product name

ZBTB33 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZBTB33 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 68-74 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006777).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116913(ZBTB33 Antibody) at dilution of 1:100



IP Result of anti-ZBTB33 (IP:Catalog No:116913, 4ug; Detection:Catalog No:116913 1:300) with A431 cells lysate 1200ug.