

ZBTB33 antibody

Catalog Number: 116913

Product name

ZBTB33 antibody

Specificity

Human; other species not tested.

Antibody description

ZBTB33 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 68-74 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006777). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:116913(ZBTB33 Antibody) at dilution of 1:100



IP Result of anti-ZBTB33 (IP:Catalog No:116913, 4ug; Detection:Catalog No:116913 1:300) with A431 cells lysate 1200ug.