

ZBTB46 antibody



Catalog Number: 116918

Product name

ZBTB46 antibody

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

ZBTB46 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of ZBTB46 recombinant protein (Accession Number: NM_001369741). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116918(ZBTB46 Antibody) at dilution of 1:600