

ZC3H12D antibody



Catalog Number: 116924

Product name

ZC3H12D antibody

Immunogen

[Human ZC3H12D Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZC3H12D Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in Raji cells, K-562 cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 55-58 kDa

Preparation

This antibody was obtained by immunization of ZC3H12D recombinant protein (Accession Number: NM_207360). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations




Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116924(ZC3H12D Antibody) at dilution of 1:600



IP Result of anti-ZC3H12D (IP:Catalog No:116924, 4ug; Detection:Catalog No:116924 1:600) with Raji cells lysate 2000ug.

ZC3H12D antibody

Catalog Number: 116924



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:116924(ZC3H12D Antibody) at dilution of 1:50 (under 40x lens)