

# ZC3H7B antibody



Catalog Number: 116927

## Product name

ZC3H7B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZC3H7B Rabbit Polyclonal antibody. Positive WB detected in HUVEC cells, mouse brain tissue. Observed molecular weight by Western-blot: 111 kDa

## Preparation

This antibody was obtained by immunization of ZC3H7B recombinant protein (Accession Number: NM\_017590). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:116927(ZC3H7B Antibody) at dilution of 1:600