

# ZC3H8 antibody



Catalog Number: 116928

## Product name

ZC3H8 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZC3H8 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, PC-3 cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of ZC3H8 recombinant protein (Accession Number: BC032001). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116928(ZC3H8 Antibody) at dilution of 1:2000