

ZC3H8 antibody

Catalog Number: 116928

Product name

ZC3H8 antibody

Specificity

Human; other species not tested.

Antibody description

ZC3H8 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, PC-3 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of ZC3H8 recombinant protein (Accession Number: BC032001). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116928(ZC3H8 Antibody) at dilution of 1:2000