

# ZNF114 antibody



Catalog Number: 116945

## Product name

ZNF114 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse; other species not tested.

## Applications

ELISA, WB

## Antibody description

ZNF114 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 40-50 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of ZNF114 recombinant protein (Accession Number: NM\_001331097). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116945(ZNF114 Antibody) at dilution of 1:600

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality