

ZNF146 antibody



Catalog Number: 116950

Product name

ZNF146 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNF146 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse lung tissue. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 20-25kd, 33-35kd

Preparation

This antibody was obtained by immunization of ZNF146 recombinant protein (Accession Number: BC005154). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116950(ZNF146 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human heart using Catalog No:116950(ZNF146 antibody) at dilution of 1:100 (under 10x lens)