

ZNF195 antibody



Catalog Number: 116957

Product name

ZNF195 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF195 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of ZNF195 recombinant protein (Accession Number: NM_001130519). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116957(ZNF195 Antibody) at dilution of 1:1000