

# ZNF2 antibody

Catalog Number: 116958

## Product name

ZNF2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZNF2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human spleen tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of ZNF2 recombinant protein (Accession Number: NM\_001017396). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116958(ZNF2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116958(ZNF2 antibody) at dilution of 1:50 (under 10x lens)