

# XBP1 antibody



Catalog Number: 116982

## Product name

XBP1 antibody

## Immunogen

[Human XBP1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

XBP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of XBP1 recombinant protein (Accession Number: NM\_005080). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116982(XBP1 Antibody) at dilution of 1:300