

XIRP2 antibody



Catalog Number: 116987

Product name

XIRP2 antibody

Specificity

human, mouse; other species not tested.

Antibody description

XIRP2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle. Positive WB detected in COLO 320 cells, human heart tissue. Observed molecular weight by Western-blot: 260-300kd

Preparation

This antibody was obtained by immunization of XIRP2 recombinant protein (Accession Number: BC022888). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116987(XIRP2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:116987(XIRP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:116987(XIRP2 antibody) at dilution of 1:50

XIRP2 antibody



Catalog Number: 116987

(under 40x lens)