

XK antibody

Catalog Number: 116988

Product name

XK antibody

WB: 1:200-1:2000

IHC: 1:50-1:400

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

XK Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, MCF-7 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 45-51 kDa

Preparation

This antibody was obtained by immunization of XK recombinant protein (Accession Number: NM_021083). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:116988(XK Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:116988(XK Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded

XK antibody



Catalog Number: 116988

human pancreas cancer tissue slide using Catalog No:116988(XK Antibody) at dilution of 1:200 (under 40x lens).