

# XLF antibody

Catalog Number: 116990

---

## Product name

XLF antibody

## Specificity

Human; other species not tested.

## Antibody description

XLF Rabbit Polyclonal antibody. Positive WB detected in human small intestine tissue. Observed molecular weight by Western-blot: 33kd

## Preparation

This antibody was obtained by immunization of XLF recombinant protein (Accession Number: NM\_001377498). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:116990(NHEJ1 antibody) at dilution of 1:300